

ALTER 

WE CREATE PRODUCTS  
THAT CHANGE THE WAY  
PEOPLE LIVE.

*And Walk. And Run.*



**ANTI-GRAVITY TREADMILL® M320**

# THE AlterG® ANTI-GRAVITY TREADMILL® M320

## DEFY GRAVITY

When your patients are recovering from surgery or injury, the AlterG® Anti-Gravity Treadmill® M320 offers a revolutionary approach to rehabilitation – the ability to defy gravity and work towards regaining function like never before. AlterG's unique differential air pressure technology provides accurate, safe and comfortable partial weight-bearing therapy while promoting normal gait patterns.

The AlterG Anti-Gravity Treadmill M320 expands the concept of “Weight Bearing as Tolerated” by offering precise unweighting in one percent increments to as low as 20 percent of the patient’s body weight. AlterG's unique precision gives patients the ability to set the precise point where exercise becomes pain free and provides clinicians a way to accurately measure patient progress. AlterG is used and recommended by leading physical therapy clinics, hospitals, colleges and pro sports teams.



### Patients Recover and Prefer AlterG

The AlterG Anti-Gravity Treadmill M320 provides precise early partial weight-bearing following injury or surgery to restore natural gait and range of motion. It enables neuromuscular and proprioceptive retraining that enhances balance and confidence, and improves gait and stride length. The AlterG also provides “prehab” conditioning with significantly reduced pain. Patients who try it overwhelmingly prefer it and 100% of those surveyed say they like it.

#### The FDA-cleared AlterG Anti-Gravity Treadmill M320 can be used for:

- Rehabilitation following injury or surgery of the lower extremity (hip, knee, ankle or foot)
- Rehabilitation after total joint replacement
- Gait training in neurologic patients
- Strengthening and conditioning in older patients
- Weight control and reduction
- Sport specific conditioning programs
- Aerobic conditioning

“This equipment will revolutionize rehab, and its arrival coincides with the increasing pressure to control the healthcare costs by significantly speeding up the recovery process.”

**Karen Brown**, PT, Owner - Alliance Physical Therapy, NE

### Athletes Train Like Never Before

The world’s best athletes and sports teams consider the AlterG Anti-Gravity Treadmill M320 an essential part of their athletic conditioning and rehabilitation programs:

- Strengthens and improves coordination of muscles, which in turn protects surrounding joints.
- Promotes the full range of motion while minimizing stress.
- Enables injured athletes to maximize their fitness retention as they recover.

### Increase Your Revenue

As the only rehabilitation and training equipment of its kind, the AlterG Anti-Gravity Treadmill M320 is redefining the practice of body weight supported therapy. Its growing recognition and popularity in the medical community can help you expand your practice through new patient referrals and additional patient visits. Many clinics have found that this new approach to training and fitness has added cash-paying patients to their clientele. Now you can lease an AlterG for a low monthly cost. When you make the choice to defy gravity, you offer your patients the chance to recover better, and train like they never have before.

MODEL	M320
<b>PERFORMANCE</b>	
- User weight range	85 - 400 lbs
- Body weight adjustment range	0 - 80% reduction*
- Body weight adjustment increment	1%
- Speed range	0 - 12 mph
- Speed adjustment increment	0.1 mph* **
- Incline range	0 - 15%
- Reverse	0-3 mph

MODEL	M320
<b>DIMENSIONS</b>	
- Length	84 inches
- Width	38 inches
- Height	74 inches
- Weight	750 lbs

\* Maximum body weight reduction is limited to 65% for user weights above 320 lbs  
 \*\* For speeds above 10 mph, increment changes to 1.0 mph